

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS16 3132

Goat anti-Dog serum proteins

Product information

Immunogen Pooled whole dog serum and partly purified serum fractions

Host Goat

Clonality Polyclonal

Purity Serum

Format Lyophilized

Quantity 1 ml

Reconstitution For reconstitution add 1 ml of sterile destilled water

Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to

the cap or sides of the tube. Reconstituted material can be stoored in 4°C up to 7 days.

Application information

Recommended dilution

In immunoelectrophoresis use 2 µl serum or equivalent against 120 µl antiserum. In double radial immunodiffusion (Ouchterlony) use a rosette arrangement with 10 µl antiserum in 3 mm diameter center well and 2 µl serum samples (neat and serially diluted in 2 mm diameter peripheral wells. Different bleedings of the immunized animals are pooled to obtain a broad spectrum balanced against the varying concentrations of the individual serum protein components.

Confirmed reactivity Dog serum proteins