

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS07 224 Fdx3 | Ferredoxin 3

Product information

Immunogen KLH-conjugated synthetic peptide chosen from a C-terminal of Chlamydomonas reinhardtii Q2HZ24

Host Rabbit

Clonality Polyclonal

Purity Immunogen affinity purified serum in PBS pH 7.4.

Format Lyophilized

Quantity 200 μg

Reconstitution For reconstitution add 166 μl of sterile water

Storage Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to

the cap or sides of the tube.

Application information

Recommended dilution 1 : 5000-1 : 10 000 (WB)

Expected | apparent 19 | 17 kDa

MW

Confirmed reactivity Chlamydomonas reinhardtii

Predicted reactivity Chlamydomonas reinhardtii

Not reactive in No confirmed exceptions from predicted reactivity are currently known

Terrauchi et al. (2009). Pattern of expression and substrate specificity of chloroplast ferredoxins from Chlamydomonas

reinhardtii. J Biol Chem 284 (38):25867-25878.